

advantages of the
intense pulsed light
microprocessed multi pulse

- ALL SKIN TONES
- ALL HAIR TYPES
- CAN BE APPLIED THROUGHOUT THE WHOLE YEAR
- ALMOST PAINLESS TECHNOLOGY
- MAXIMUM SAFETY AND EFFECTIVENESS



+ PHOTOEPILATION
FACIAL PHOTOREJUVENATION



termosalud.com



© 2014 Termosalud - Marketing Dept.



The pulsed light
that protects your skin



+ PHOTOEPILATION
FACIAL
PHOTOREJUVENATION





What is AURA?

AURA is an intense pulsed light (IPL) equipment developed to carry out photoepilation and facial photorejuvenation treatments.

The IPL system manufactured by Termosalud, is capable of generating a multitude of sub-pulses with each shot of light, which makes the treatment more effective and gives an excellent level of safety.


What is the intense pulsed light, IPL?

Intense pulsed light is an advanced, non-invasive and painless technology which makes possible to treat two very common aesthetic problems.

It is a very effective technique for the **removal of unwanted hair**; it is also used very successfully for the treatment for the **skin ageing** and benign anti-aesthetic conditions, such as blemishes caused by age or the sun and congenital vascular and pigmentation conditions.

Its **“all types” technology** generates and controls the energy in a completely customised way. This technology is effective on **all skin tones, all hair types and all year round**.


It is a versatile system suited for every person and every problem, which allows a **combined photoepilation and facial rejuvenation treatment**.



PHOTOEPILATION IPL

IPL photoepilation is a **permanent and long-lasting** hair removal system, suited for all skin types, using a flash of light emitted by a high power Xenon lamp using **Intense Pulsed Light (IPL)** equipment.

The melanine absorbs the light and carries it to the follicle, where it is transformed into heat, burning the germinative cells and causing **the pilous bulb to be destroyed** and therefore destroying its capacity to generate a new hair, all without damaging the surrounding tissue at all.



After the treatment, the hair will fall out naturally through a physiological process in approximately three weeks.

TERMO SALUD




FACIAL PHOTOREJUVENATION

As we get older the skin becomes thinner and drier. The lower collagen makes the skin less elastic and the first wrinkles start to appear.

IPL photorejuvenation consists of applying pulses of intense light, with a defined range of wavelengths which penetrate the skin, heating it deep under the surface and stimulating the fibroblasts causing them to produce new collagen and elastin, thereby **improving the quality of the skin** and giving it a long-lasting smooth, firm look.

REGENERATION OF THE SKIN:



Another benefit of this phototherapy is that it **evens out the tone of the skin**, reducing the chromatic differences that can be caused by melanin (sun spots) or can be **vascular**.